

# HS60 for CV

Focus on your needs

## Focus on extended diagnostic capabilities

Extend diagnostic boundaries with a versatile ultrasound HS60. It enables users to diagnose with accuracy various anatomies including cardiac and vascular using dedicated app and advanced imaging features.



### Product information

- S-Vision™ imaging engine
- S-View™ transducer technology (PA1-5A)
- Advanced imaging functions (S-Harmonic™, ClearVision)
- Cardiac solutions (Strain+, Stress Echo)
- Semi-automated measurement of intima-media thickness (Auto IMT+)
- Workflow improving tool (Quick Preset, EZ-Exam+™, EZ-Compare™, Advanced QuickScan)
- 21.5-inch LED monitor / 10.1-inch touch screen / Solid State Drive / Gel warmer



## Comprehensive selection of transducers

### Linear array transducers



LA4-18BD



LA3-16A



LA3-14AD



LA2-9A

### Phased array transducers



PA1-5A



PA3-8B


### CW transducers




DP2B




CW6.0



Strain+



Auto IMT



Mitral valve with PW



4 chamber view